




P.O. Box 222, Lebanon, OH 45036 Phone: 513.932.7379 Fax: 513.932.4493


Approximate dimensions of STANDARD Tri-Clamp style gaskets (40MP-)

DIMENSIONS	ID	BEAD	OD	
1/2"	.390		.855	For 1/2" and 3/4" Approx. thickness: .079" Approx. bead height: .180"
3/4"	.645		.855	
1"	.870	1.718	1.975	For 1" through 12" Approx. thickness: .065" (OD) and .085" (ID)
1 1/2"	1.385	1.718	1.975	
2"	1.885	2.218	2.505	
2 1/2"	2.371	2.781	3.035	
3"	2.865	3.281	3.570	
4"	3.855	4.345	4.670	
5"	4.795	5.281	5.670	
6"	5.795	6.180	6.555	
8"	7.795	8.156	8.555	
10"	9.705	10.156	10.545	
12"	11.735	12.130	12.535	



Approximate dimensions of FLANGED Tri-Clamp style gaskets (40MPF-)

DIMENSIONS	ID	BEAD	OD	
1/2"	.390		.855	For 1/2" and 3/4" Approx. thickness: .079" Approx. bead height: .180"
3/4"	.645		.855	
1"	.870	1.718	1.200	For 1" through 12" Approx. thickness: .065" (OD) and .085" (ID)
1 1/2"	1.385	1.718	1.200	
2"	1.885	2.218	2.530	
2 1/2"	2.371	2.781	3.060	
3"	2.865	3.281	3.595	
4"	3.855	4.345	4.695	
5"	4.795	5.281	5.695	
6"	5.795	6.180	6.580	
8"	7.795	8.156	8.580	
10"	9.705	10.156	10.570	
12"	11.735	12.130	12.560	



Material	FDA Certified		USP Class VI Certified	
	(Standard TriClamp)	(Flanged TriClamp)	(Standard TriClamp)	(Flanged TriClamp)
Buna-N	40MP-U	40MPF-U	(not available in Class VI)	
EPDM	40MP-E	40MPF-E	40RX-E	40RXF-E
Viton	40MP-SFY	40MPF-SFY	40RX-SFY	40RXF-SFY
Silicone	40MP-XPL	40MPF-XPL	40RX-PX	40RXF-PX
Teflon	40MP-G	N/A	40RX-G	N/A
Silverback*	40MP-GST	N/A	40RX-GST	N/A

Silverback material is a 50/50 blend of stainless steel fibers and PTFE